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Storage information

The storage location should be dry and at room temperature. The product has to be protected against direct sunlight exposure and humidity.

Cleaning information

Clean the brace with moist cloth using mild detergents.

Material

Textile: Polyester, elastane, cotton, silicone, elastomer.

Guarantee

The guarantee lies within the framework of the legal regulations.

Disposal

The product can be disposed of with the normal household waste.

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 mdh sp. z o.o.
ul. Maratońska 104
94-007 Łódź
www.qmedbymeira.pl
ortopedia@mdh.pl
tel. +48 42 636 81 44

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EN

As-15 Medical arm sling with abduction cushion

Thank you very much for selecting a product of the Meyra-Ortopedia Group. The product complies with all legal regulations according to the Medical Products Law.

Please read the directions for use carefully and keep to the advice and recommendations given. Should you have further questions please address the responsible physician, your local specialist dealer or our service hotline.

Product description

The abduction cushion with arm sling is used for the treatment of different indications in the area of shoulder, upper arm and elbow. The fore arm steadily rests in a padded arm ling which is fixed around shoulder and waist with quick-lock buckles. For a higher convenience the shoulder strap is padded. The abduction cushion spreads the arm from the body as required after e.g. rotator cuff surgery. This position hastens the healing process.

Safety instructions

The first application of the product must be carried out by special personnel trained on these kind of products as other-wise the product liability cannot be guaranteed. Basically, when adjusting the Velcro straps it has to be ensured that they are not fastened too tight to avoid constrictions. The product must not be worn on open wounds.

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Indications

The abduction cushion is used after rotator cuff surgery, injuries of ligaments, tendons and muscles on shoulder, upperarm and elbow.

Application

1. Adjust the shoulder strap to wide length and close it.
2. Lay the forearm into the arm rest, taking the thumb into the body side loop and the fingers into the loop on the body far side.
3. Place the fixation straps over head and side of the healthy shoulder and adjust them for a good fit without constrictions.
4. The forearm strap is laid around the waist of the armrest and adjusted to the required length.
5. The semicircular side of the abduction cushion is positioned to the body side and the outer Velcro strap is connected with the Velcro strap of the armrest.

Contraindications

Provided that the product is employed properly contraindications are not known. If the slings are too tight pressure symptoms may evolve.

Employment according to purpose and indications

The abduction cushion is used in the case of treatment after dislocation, arthroscopy, and for general stabilization of the shoulder joint.

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Further information

- The product should only be worn under the care of a physician and for the intended clinical pictures.
- The product is intended for the treatment of one patient.
- Improper application excludes the product liability of the manufacturer.
- Should the symptoms aggravate a physician should be consulted immediately.
- A combination with other products has to be agreed with the physician in advance.

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