PRODUCT DESCRIPTION:

Spincare Flex II is made of flexible straps. The fabric carries humidity and warmth away. Additional knitwear applied in the dorsal part improves comfort of use and two oblique bands support dorsal muscles. System of additional straps enables easy adjustment of dorsal compression and control of dynamic posture correction. Four semi-elastic, asymmetric stays provide passive stabilization and enable optimal body adjustment.

PUTTING ON:

- Do not forget to make sure that the brace size has been chosen carefully.
- In order to put on the brace, the patient should lie down.
- Put the brace so that it would be placed 2 cm above the coccygeal vertebrae - the brace's label (you can find it inside the brace) should be directed upward.
- The back lines should adhere to T-L-S area.
- Fasten the belts and regulate the dorsal pressure.
- Make sure that the brace provides optimal stabilization but does not disturb breathing.
- The time of using the brace should be precised by medical specialist and should be accompanied by rehabilitation exercises program.

| Size | 1 | 2 | 3 | 4 | 5 | 6 |
|--------------------------|---------|---------|---------|----------|-----------|-----------|
| Waist circumference (cm) | 52 - 62 | 63 - 74 | 75 - 89 | 90 - 105 | 106 - 120 | 121 - 140 |

CONSERVATION AND MAINTENANCE:

The brace should be washed with gentle detergents in max. temperature 35°. Do not chlorinate, do not dry clean, wring gently before drying. Do not dry close to the artificial heat sources. Dry in horizontal position. Do not iron.

ATTENTION!

Not observing the rules of personal hygiene and cleaning the vest with inappropriate detergents may result in skin abrasion or other effects related to perspiration. Please, do not use gels or warming ointment. Always rinse carefully after washing-detergent residue can cause skin irritation and damage product material.
The product is not fireproof. Do not use close to the open sources of fire.

THE FINAL DECISION CONCERNING APPLICATION AND CHOICE OF THE EQUIPMENT BELONGS TO A MEDICAL SPECIALIST.

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INSTRUCTIONS FOR USE

SPINCARE FLEX II – thoracic-lumbosacral orthosis

Spincare Flex II brace helps effectively with rehabilitation and treatment in case of injuries, diseases and dysfunctions of intervertebral joints, intervertebral disc, muscle attachments, ligaments of lumbar area and ligaments of lumbosacral passage. The dynamic stabilization enables effective control of movement range, relief and posture correction which alleviate symptoms connected to degenerative changes and functional disorders of spinal segment.

INDICATIONS:

- Non-operative treatment of spinal segment dysfunctions of different etiology.
- Stabilization of lumbar area in case of segmental instability after surgeries.
- Follow-up treatment after neurosurgical/orthopedic interventions.
- Prophylaxis against injuries, overload.
- Congenital defects of lumbar area (with or without pain syndromes).
- Degenerative spinal disorders with pain syndrome.