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PLEASE READ CAREFULLY USING
BEFORE ALL INSTRUCTIONS

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
Introduction

Thank you for choosing our product. We are strongly believe it will meet your expectations. Please get acquainted with this document before first use. This instruction is delivered to give necessary information about safe use of the device. Obeying the recommendations will influence the durability and aesthetics of the product. We would like to indicate that we keep trying to improve the quality of our products which may cause minor changes that may not be included in the instruction.

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1. Knob regulating the resistance of the pedals
2. Rehabilitation pedals
3. Bicycle shoulders
4. The main frame
5. Lock of submission
6. Base arms

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Application

Rehabilitation and recreation, recommended for upper and lower limbs paralysis and orthopaedic diseases.

Exercises with use of the exerciser allow to form the muscle power, improve blood circulation within limbs, stimulate nervous system and improve the joint range of motion. Especially recommended for exercises during paralyse of lower limb peripheral nerves and for those who have experienced the cerebral stroke or other neurological dysfunction.

The possibility of applying the exerciser at hospitals as well as in household conditions.

Contraindication

Conditions requiring limited upper and lower limbs motion.

contraindication to physical movements (high blood pressure, cardio-vascular and breathing disorder)

Use

UNFOLDING

- After removing the packaging, make sure the exerciser is not mechanically damaged.
- Then straighten the base frame until the lock of submission is latched in the main frame.
- Make sure the legs are properly installed and that do not move.
- Knob governing the resistance of the pedals lightly screw into the hole located in the central part of the frame.

ASSEMBLING

- To submitted the exerciser, press lock of submission and submit his arms toward the main frame.

Safety of use WARNING!

- a) Before using exerciser for upper and lower limbs , refer to this manual and keep it for future reference.

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Everyday use and maintenance

LOWER AND UPPER LIMBS EXERCISE

- Sitting on a chair, set the exerciser in front of you, with the pedals located within the user's legs.
- Sitting on a chair, set the exerciser in front of you on a table or other piece of furniture with a flat area with that the pedals located within the user's hand.
- Next, insert your feet / hands into pedals stirrups. if you feel discomfort caused by improper distance between the pedals and the chair, it is recommended to set a proper distance.
- Now you can start exercises.

ADJUSTMENT

- In the central part of the exerciser is a knob which regulate the resistance of the pedals.
- When You turn the knob clockwise, the resistance is increased , it means that you should use more force to move the pedals.
- When You turn the knob counterclockwise, the resistance is decreased , it means that you should use less force to move the pedals.

CALORIE COUNTER INSTALLATION ONLY FOR MODELS WITH CALORIE COUNTER

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- Before installing the calorie counter, unscrew the knob regulating the resistance of the pedal.

- Uchwyt mocowania
Montážní držák
- Naklejka informacyjna
Informační nálepka
- Zaczep
Háček

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- The calorie counter should be installed in such a way that information sticker was on the same side like of the sensor of counter.
- Thread must be surrounded by the mounting bracket . Knob regulating the resistance of the pedals must be screwed in the thread. Hooks should cover the main frame of exerciser. (Illustration 4.)

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BATTERY REPLACEMENT

- The calorie counter is operated by one battery AAA 1,5V.
- When the battery is at the finish must be replaced.

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EVERYDAY USE AND MAINTENANCE

- Clean the main frame, arms, pedal shoulders with a wet cloth.
- Systematically remove the dirt from the movable elements of the rotor (screws, resistance mechanisms).

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You must not use:

- solvents,
- toilet detergents,
- sharp brushes and hard objects,
- detergents that contain chlorine,
- corrosive detergents,
- Do not clean the exerciser with pressure or steam aggregates!
- For polished and chrome parts care it is recommended to use detergents for polished objects that are available on the market.

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Frequency	Object of control	Comments
Before first installation	Checking the technical condition and all the elements within the packaging.	Check if any elements have any visible mechanical damage.
Every two weeks	Screw and nut control regarding their correct fixing.	Make sure if all the screws and nuts are fixed tightly.
Every two months	Cleanliness and general condition.	In case of mechanical damage immediately contact the producer.

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Repair

- In case of reparation needs contact the producer.

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
Customer service

- Should you have any questions or need help, contact the producer.

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