ROM ADJUSTMENT:

The changes and limiting of angular movement is prescribed by a doctor. If your doctor tells you to limit the angular range of motion: Loosen two screws on the dial. Loosen the two screws on the dial if you want to change the range of motion. Remove the wedge held by the screw. Replace it with a wedge with appropriate angular range marking. When the range of flexion and extension angle is set within the therapeutic range, tighten the screw. Note that the angular range should be determined before the orthosis is applied, and that both dials have to be set to the same value.

Size	XS	S	M	L	XL	XXL	XXXL**
Hip circumference (cm)	30 - 36	36 - 45	44 - 50	50 - 57	57 - 64	64 - 72	72 - 78

* The size is selected based on the circumference of the thigh (in cm), measured 15cm above the centre of the patella

** open version only

MAINTENANCE:

Remove the liners before washing. Hand wash in warm water with mild detergent. Do not use chlorine. Do not wring, gently squeeze. Dry in the extended position away from any artificial heat sources. Do not iron.

ATTENTION!

Not observing the rules of personal hygiene and cleaning the vest with inappropriate detergents may result in skin abrasion or other effects related to perspiration. Please, do not use gels or warming ointment. Always rinse carefully after washing- detergent residue can cause skin irritation and damage product material.

The product is not fireproof. Do not use close to the open sources of fire.

THE FINAL DECISION CONCERNING THE APPLICATION AND CHOICE OF THE EQUIPMENT BELONGS TO A CONSULTANT.

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REFLEX LONG – Knee Brace

INDICATIONS:

- Mild instability of the knee joint of various etiologies,
- Conservative treatment of knee ligaments injuries and sprains,
- Degenerative joint disease or osteoarthritis of the knee joint.

CONTRAINDICATIONS:

- Conditions requiring immobilization, long- healing wounds and ulcers, significant bone deformation.

ATTENTION! The brace should be used according to the instructions given by the doctor or physiotherapist. When used for the first time, the support should be put on by the experienced health-care personnel. Range of motion should be adjusted according to instructions given by doctor or physiotherapist. Before the first use, please, read the user's manual carefully. If you feel discomfort while using the brace, in case of swelling or other skin problems, please consult a doctor or physiotherapist. To increase the efficiency of the orthosis, adjust the shape of the stays and the position and length of the knee stabilizing ring to individual need. Before fitting the orthosis, remove the stays from the pockets in the body. Adjust stays so that they fit the limb shape. Place the rails back in the pockets and make sure they are properly fastened.

FITTING:

- Closed version: Pull the brace through the foot with the help of belt loops. Kneecap should be in a dedicated hole.
- Open version: Fasten the brace around the knee starting from the lower part. Kneecap should be in a dedicated hole.