CONSERVATION AND MAINTENANCE:

The brace should be washed using gentle detergents in max. temperature 35°. Do not chlorinate, do not dry clean, wring gently before drying. Do not dry close to the artificial heat sources. Dry in horizontal position. Do not iron.

ATTENTION!

Not observing the rules of personal hygiene and cleaning the vest with inappropriate detergents may result in skin abrasion or other effects related to perspiration. Please, do not use gels or warming ointment. Always rinse carefully after washing-detergent residue can cause skin irritation and damage product material.

The product is not fireproof. Do not use close to the open sources of fire.

THE FINAL DECISION CONCERNING THE APPLICATION AND CHOICE OF THE EQUIPMENT BELONGS TO A MEDICAL SPECIALIST.

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INSTRUCTIONS FOR USE

SPINCARE EXTENSOR - Upright keeper

Spincare Extensor is an upright keeper of wide therapeutic use. It effectively helps rehabilitation and treatment in case of collar bone injuries and facilitates correction of postural defects. It can be also used as a preventive aid to establish the habit of correct body posture.

INDICATIONS:

- Collar-bone injury which does not require surgical or orthopaedic intervention.
- Posture defects which do not require surgical intervention.
- Pain syndromes of thoracic section.
- Kyphosis.

PRODUCT DESCRIPTION:

Spincare Extensor is made of flexible, tensile straps. The fabric carries humidity and warmth away. Four semi-elastic stays reinforce dorsal part. Reinforcing pads placed on the shoulder straps provide proper adhesion and prevent them from creasing which is crucial for user's comfort. Special construction and wide range of sizes allow to correct the body posture in patients at different ages.

PUTTING ON:

1. Unfasten the Spincare Extensor and put it on the thoracic area so that the dorsal part would be placed between the scapular bones.
2. Put the straps on the shoulders, under the armpits and cross them at the back, Then lace them through belt loop of the dorsal pad.
3. Fasten the straps on the chest, below the sternum edge.

| Size | 1 | 2 | 3 | 4 | 5 | 6 |
|--------------------------|---------|---------|---------|---------|---------|----------|
| Chest circumference (cm) | 40 – 50 | 50 – 60 | 60 – 70 | 70 – 80 | 80 – 90 | 90 - 100 |