Qmed®
Active

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Available in four resistance levels.

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Hmotnost: 80g

EN

The mass is ideal for the rehabilitation of hands and fingers. Regular exercise reduces pain and stiffness in the joints of the fingers, hand, and wrist, as well as improves range of motion in the interphalangeal joints and wrist, and helps to strengthen your grip or restore impaired gripping functions. It is perfect for rehabilitation in the case of degenerative diseases of the hand, after injuries, with weakened strength of the forearm and hand muscles.

Hand exercises using the Therapy Putty mass consist in rolling the material with the outer and inner sides of the hand, crushing, squeezing or pressing the mass with single fingers. The mass is made of non-toxic silicone. Keep away from textiles.

Available in four levels of hardness, differing in color:

- Extra light (green)
- Light (yellow)
- Medium (red)
- Heavy (blue)

Recommendations:

- Rheumatoid arthritis, arthritis, carpal tunnel syndrome, tendinitis, tennis elbow, recovery from finger or wrist fracture or tendon surgery.
- Complementing the training session with exercises using hand and finger weight will help you avoid injuries and reduce the muscle imbalance (imbalance) caused by excessive hand and finger exercises.
- Conditions requiring increased finger and grip strength

Limitations:

- Do not use in patients with impaired sensation and people with skin lesions such as ulcers, blisters, abrasions, etc.
- The occurrence of pain during exercise
- Conditions immediately after sprains or bone fractures
- Incomplete bone union
- Inflammation of joints and periarticular tissues
- Swelling
- Fresh postoperative scars or other wounds

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Make a hook fist and squeeze putty between your fingers and your hand.

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
Form a ball with the putty. Squeeze the putty with all four fingers.

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Loop the putty over the end of the finger while it's bent. Straighten your finger using the resistance of the putty.

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
Spread the putty like a pancake over your fingers and thumb. Try to spread them apart.

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Place a 1 inch thick piece of putty between each pair of fingers and squeeze together.

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Pinch the putty between each finger and the thumb.

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Put a tight ring of putty around two fingers. Spread your fingers apart. Repeat the exercise for each finger.

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
Loop the putty at the end of your thumb while it is bent. Try to straighten your thumb by pulling it upward.

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Push your thumb into the putty as you move the thumb toward your small finger.

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Press the putty with your thumb against the side of your index finger. Keep your fingers and thumb straight.

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Squeeze the putty between your thumb and side of your index finger.

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Pull the putty using your thumb, index and middle fingers.